

Ocean Optics SeaBreeze UML


This presentation provides an overview of the object and event models for the Ocean Optics SeaBreeze device driver. This information is Confidential/Proprietary Ocean Optics, Inc. Further redistribution is prohibited.


Copyright © 2009 Ocean Optics, Inc. All Rights Reserved.


Note that, due to space constraints, not all classes are shown in these UML diagrams; a representative subset may be given where the entire set would not easily fit. Source code should be provided along with this document, and that can be used to determine the extent of the classes included in SeaBreeze.













- + readSpectrum (Bus bus) : vector<double> *
- + readUnformattedSpectrum (Bus bus): vector
byte> *
- + requestSpectrum (Bus bus) : void
- + setIntegrationTimeMicros (Bus bus, unsigned long t): void

OOISpectrometerProtocol

- # integrationTimeExchange : ooiProtocol::IntegrationTimeExchange *
- # requestSpectrumExchange : Transfer *
 # spectrumTransferExchange : Transfer *
- # unformattedSpectrumExchange : Transfer *
- + readSpectrum (Bus bus) : vector<double> *
- + readUnformattedSpectrum (Bus bus) : vector
byte> *
- + requestSpectrum (Bus bus) : void
- + setIntegrationTimeMicros (Bus bus, unsigned long t): void

unformattedSpectrumExchange ↓ 1

spectrumTransferExchange 1

ooiProtocol::ReadSpectrumExchange

. . .

ooiProtocol::IntegrationTimeExchange


...


request Spectrum Exchange


1


request Spectrum Exchange

ooiProtocol::RequestSpectrumExchange









Feature

protocols : vector<Protocol *>

lookupProtocolImpl (Protocol) : ProtocolHelper *

OOISpectrometerFeature

integrationTimeBase : long
integrationTimeIncrement : long
integrationTimeMaximum : long
integrationTimeMinimum : long

numberOfPixels : int
maxIntensity : int

- + getIntegrationTimeIncrement (): long
- + getIntegrationTimeMaximum (): long
- + getIntegrationTimeMinimum (): long
- + getMaximumIntensity (): int
- + getNumberOfPixels () : int
- + getSpectrum (Protocol proto, Bus bus) : vector<double> *
- + getUnformattedSpectrum (Protocol proto, Bus bus) : vector
byte> *
- + getWavelengths (Protocol proto, Bus bus) : vector<double> *
- + readUnformattedSpectrum (Protocol proto, Bus bus): vector
byte> *
- + setIntegrationTimeMicros (Protocol proto, Bus bus, unsigned long t) : void
- + writeRequestSpectrum (Protocol proto, Bus bus) : void

QE65000SpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION_TIME_INCREMENT : long
- INTEGRATION_TIME_MAXIMUM : long
- INTEGRATION_TIME_MINIMUM : long

HR2000PlusSpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION_TIME_INCREMENT : long
- INTEGRATION_TIME_MAXIMUM : long
- INTEGRATION_TIME_MINIMUM : long

MayaPro2000SpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION_TIME_INCREMENT : long
- INTEGRATION_TIME_MAXIMUM : long
- INTEGRATION_TIME_MINIMUM : long

HR2000SpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION_TIME_INCREMENT : long
- INTEGRATION_TIME_MAXIMUM : long
- INTEGRATION_TIME_MINIMUM : long

Maya2000SpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION TIME INCREMENT: long
- INTEGRATION_TIME_MAXIMUM : long
- INTEGRATION TIME MINIMUM: long

HR4000SpectrometerFeature

- INTEGRATION_TIME_BASE : long
- INTEGRATION_TIME_INCREMENT: long
- INTEGRATION TIME MAXIMUM: long
- INTEGRATION_TIME_MINIMUM : long

NIRQuestSpectrometerFeature

- # INTEGRATION_TIME_BASE : long
- # INTEGRATION_TIME_INCREMENT : long
- # INTEGRATION_TIME_MAXIMUM : long
- # INTEGRATION_TIME_MINIMUM : long

NIRQuest512SpectrometerFeature

NIRQuest256SpectrometerFeature



protocols : vector<Protocol *>

lookupProtocolImpl (Protocol) : Protocol *

EEPROMSlotFeatureBase

readEEPROMSlot (Protocol proto, Bus bus, unsigned int slot) : vector<byte> *
readDouble (Protocol proto, Bus bus, unsigned int slot) : double

EEPROMSlotFeature

- numberOfSlots : unsigned int
- + readAllEEPROMSlots (Protocol proto, Bus bus) : vector< vector

 byte> *> *
- + readEEPROMSlot (Protocol proto, Bus bus, unsigned int slot): vector
byte> *

SerialNumberEEPROMSlotFeature

+ readSerialNumber (Protocol proto, Bus bus) : string *

WavelengthEEPROMSlotFeature

- numberOfPixels : unsigned int
- + readWavelengths (Protocol proto, Bus bus) : vector<double> *
- # computeWavelengths (double polynomial[], int length) : vector<double> *

Feature

protocols : vector<Protocol *>

lookupProtocolImpl (Protocol) : Protocol *


ThermoElectricQEFeature

- + getDefaultSetPointCelsius (Protocol proto, Bus bus) : double
- + getDefaultTECEnable (Protocol proto, Bus bus) : bool
- + getTemperatureCelsius (Protocol proto, Bus bus) : double
- + setTECEnable (Protocol proto, Bus bus, bool enable) : void
- + setTemperatureSetPointCelsius (Protocol proto, Bus bus, double degC) : void
- readTECDefaults (Protocol proto, Bus bus) : vector


 byte> *

StrobeLampFeature


+ setStrobeEnable (Protocol proto, Bus bus, bool enable) : void




HR2000 Spectrum Request (simple write)



HR2000 Spectrum Read (simple read)



HR2000 EEPROM Read (write/read)



HR2000 Set Lamp Enable (value write)

